



Guía de broca dinámica combinada 2,5 mm / 3,5 mm

Combined dynamic drill guide 2.5 mm / 3.5 mm

La guía de broca dinámica brinda versatilidad para realizar orificios excéntricos compresivos y neutros para tornillos de 3,5 mm.

The dynamic drill guide provides versatility to drill eccentric compression holes and neutral holes for 3.5 mm screws.



www.saimedicine.com

[f](#) [@](#) [in](#) [v](#) [t](#) @saimedicineok

Guía de broca dinámica combinada

2.5 mm / 3.5 mm

Combined dynamic drill guide 2.5 mm / 3.5 mm





www.saimedicine.com