



Arpón Sub Knotless®

Anclaje impactado sin nudos

Knotless impacted anchor

El Sistema SUB KNOTLESS permite una reparación artroscópica simple y más segura en la inestabilidad glenohumeral.


Esta técnica Knotless elimina las complicaciones causadas por los nudos diferenciándose de otros anclajes convencionales permitiendo al cirujano colocar la sutura de forma independiente.

La tensión del tejido puede visualizarse y ajustarse, si fuera necesario, antes de la colocación final del ancla.

The Sub Knotless system allows simple and safer arthroscopic repair in glenohumeral repair.

This knotless technique eliminates the complications caused by knots, differing from other conventional anchors, allowing the surgeon to use the suture independently.

The tension of the tissue can be visualized and adjusted, if necessary, before final anchor placement.



Arpón Sub Knotless®

Sub Knotless® Anchor

Características Principales

Main Features

- + Colocación por impacto
Impact placement.
- + Alta resistencia.
High resistance.
- + Máxima resistencia a la tracción.
Maximum tensile strength.
- + Fabricado en Peek.
Made in peek.
- + Sutura ultrarresistente.
Ultra-resistant suture.
- + Iniciador descartable.
Disposable initiator.
- + Producto estéril.
Sterile product.

Indicaciones

Indications

Mano y muñeca

Hand and wrist

- + Complejo de fibrocartilago triangular
Triangular fibrocartilage complex (TFCC)
- + Ligamento colateral cubital
Ulnar collateral ligament (UCL)
- + Ligamento del escafolunar
Scapholunate ligament (SL)
- + **Hombro:** Slap / Inestabilidad
Shoulder: Shoulder Slap / Instability
- + **Cadera:** Lesión de Labrum
Hip: Hip Labrum lesions
- + **Rodilla:** Ligamento Patelofemoral
Knee: Knee Patellofemoral ligament
- + **Tobillo:** Reparación Ligamentaria
Ankle: Ankle Ligament repairing
- + **Codo:** Ligamento colateral medial / Biceps
Elbow: Elbow Collateral medial ligament / Biceps



Arpón Sub Knotless®

Sub Knotless® Anchor

3351-000

Arpón Sub Knotless 2.5

Sub Knotless Anchor Screw 2.5 mm

1901-000

Arpón Sub Knotless 2.8

Sub Knotless Anchor Screw 2.8 mm

1902-000

Arpón Sub Knotless 3.5

Sub Knotless Anchor Screw 3.5 mm

1903-000

Arpón Sub Knotless 4.5

Sub Knotless Anchor Screw 4.5 mm





www.saimedicine.com